

8/11/08 PCR V. VIMBINA #11

CDNA - 1

10x Buff - 5 x 2.5 = 12.5

(10mM) dNTPs - 1 x 2.5 = 2.5

P<sub>F</sub> - 0.5 x 2.5 = 1.25

P<sub>R</sub> - 0.5 x 2.5 = 1.25

AMPLITUDE - 0.4 x 2.5 = 1

H<sub>2</sub>O - 41.6 = 2.5 = 104

50

ADD 49 μl TO EACH TUBE

PRIMERS: SP4-SCR-Pro F/R

QPY-SP13-F/R

SCR-Pro B F/R

qPCR - VIMBINA/OYSTER SEEDINGS (16s)

CDNA - 0.5

2x IMMOMIX - 12.5 x 8 = 100

P<sub>F</sub> - 0.1 x 8 = 0.8

P<sub>R</sub> - 0.1 x 8 = 0.8

SYTO13 - 1 x 8 = 8

H<sub>2</sub>O - 10.3 x 8 = 82.4

25 μl

ADD 24 μl TO EACH WELL

8/19/08 - KLENOW OF PCR PRODUCTS FROM 8/6/08

DNA - 3

BOXERS #2 - 2.5 x 6.5 = 16.25

100μM dNTPs - 8.25 x 6.5 = 53.4

KLENOW - 0.5 x 6.5 = 3.25

H<sub>2</sub>O - 10.75 x 6.5 = 4.88 69.9

25 μl

ADD 22 μl TO EACH TUBE

INCUBATE 15min @ 25°C

ADD EDTA TO C<sub>f</sub> = 10mM

INCUBATE 75°C 20min